

The Invacare Jasmine Mobile Lift is a flexible solution that offers high comfort for all clients. The lift covers all needs, providing safe patient handling with up to 200 kg weight capacity. With its expansive, wide opening base and high lifting mechanism, Jasmine makes life easy for caregivers operating this Invacare lift system.



# **Jasmine**

The Invacare Jasmine is exceptionally mobile and easy to move around. The wide-opening base allows easy access for a wheelchair or armchair, whilst the high lifting mechanism is essential for helping clients on and off high toilet chairs or examination couches. Flexible and with optimized maneuverability, Jasmine covers all needs, at all times – for both the caregiver and the client.

Jasmine secures the patient's need for a safe and comfortable lift with control box with diagnostic system: an intelligent service monitoring system that highlights the battery status and service condition of the actuator to ensure safe patient handling. The control box provides a safe and comfortable lift for the client via 'soft start' and 'soft stop' button mechanisms.

### **Features**



Full lift range and good clearance for the user.



Electric leg adjustment. Easy and effortless solution for positioning lifter when transferring a user.



Intelligent monitoring system.



2 point spreader bar as standard.

## **Tech Specs**



Length 134cm



Min. base widt 69cm



Max. base width



Arm height 54 - 192cm



Ground clearance



Max. user Weight 200kg



Unit weight 43kg



Battery time 40 cycles

#### Warranty:

- 3 years frame
- 2 years electronics and motors
- 1 year battery



### FOR FURTHER INFORMATION:

New Zealand: Freephone 0800 468 222, email sales@invacare.co.nz or visit www.invacare.co.nz Australia: Freephone 1800 460 460, email sales@invacare.com.au or visit www.invacare.com.au Asia: Phone +61 2 8839 5330, email exportsales@invacare.com.au or visit www.invacare.com Thailand: Phone +66 0 2108 8555 or visit www.invacare.com

